

URINE FOR RESEARCH OF HAEMATOBIUM SCHISTOSOMA (bilharziasis)

MATERIAL NEEDED FOR COLLECTION

No. 3 sterile, wide-opening, screw-top containers with suction probe provided by the Collection Center or available for purchase at the pharmacy.



DIRECTIONS AND METHOD OF COLLECTION

Specimen collection must take place at the Collection Center between 11:30 am and 12:30 pm, without having urinated for at least 3 hours.

Before sample collection, the patient will be asked to climb a flight of stairs or do some knee bending or jumping.

He/she should then collect urine in the three jars by dividing the urination as follows:

1st jar: collection of the first micturition

2nd jar: collection of the middle micturition

3rd jar: collection of the last part of the micturition.

In total, about 50 mL of urine should be collected.